

Certified compliance ... The proof

Save time and buy in complete confidence by selecting furniture conforming to current regulations and standards.

Supported by tests, the **NF Santé** mark (furniture for medical/hygiene use) provides you with objective proof that the product complies with specified requirements and takes into account the appropriate criteria for this specific use in terms of:

- > **Safety**
- > **Robustness**
- > **Durability**
- > **Hygiene**

For you, it is also the assurance that the manufacturer knows how to ensure a constant level of quality for all furniture marked **NF Santé**.

WHEELCHAIRS, SEATS


- > **Safety**
 - Rearward, side and forward stability of the seat and footrests
 - Effectiveness of the brake on a 5° inclined plane
 - Fire resistance
- > **Comfort**
 - Suitable dimensions
 - Bilateral and accessible controls
 - Rolling forces
- > **Strength**
 - Thrust of 700 N on the seat back
 - Load on the seat base of 270 kg
 - Thrust of 750 N on the brake
 - Thrust of 1500 N on the footrest
- > **Durability**
 - 75 kg weight dropped 40,000 times onto the seat base
 - Seat back operated 20,000 times
 - 30,000 operations of a height-adjustable seat base when loaded with 30,000 thrusts of 400 N on the seat back
 - 30,000 operations of the footrest
 - 10,000 rolling operations of a loaded wheelchair
 - 10,000 operations of the brake
 - 40,000 alternating thrusts of 900 N on the armrests

CHEST OF DRAWERS

- > **Safety**
 - Stability with opening of loaded drawers
- > **Strength**
 - Drawers loaded to 0.3 kg/dm³
 - Load on the top of 1.5 kg/dm²
 - 5 drops from a height of 100 mm
- > **Durability**
 - Drawers opened 50,000 times
 - 2000 lateral thrusts of 300 N

CABINETS

- > **Dimensions**
 - Stability when doors opened
- > **Safety**
 - Load on pull-out shelves and on top
 - 5 drops from a height of 100 mm
- > **Durability**
 - Doors opened 40,000 times
 - 2000 lateral thrusts of 300 N testing rigidity



BEDHEAD

- > **Safety**
 - Stability on opening doors, drawers and pull-out shelves
- > **Strength**
 - Load on the drawers
 - Load on pull-out shelves and on top
 - 5 drops from a height of 100 mm
- > **Durability**
 - Drawers opened 40,000 times
 - 2000 lateral thrusts testing rigidity

BEDROOM TABLE

- > **Safety**
 - 80 kg stability for a 4-leg table
 - 35 kg for an occasional table
- > **Strength**
 - Load on the top of 1.5 kg/dm²
 - Pressure on floor less than 3 MPa
 - 5 drops by tipping onto front
- > **Durability**
 - 10,000 lateral thrusts of 300 N at the 4 corners for a 4-leg table
 - 5000 torsional thrusts of 200 N for an occasional table

TABLE FOR EATING IN BED

- > **Dimensions**
 - Stability when loaded by 300 N
- > **Strength**
 - 5 drops by tipping onto front
- > **Durability**
 - 2000 operations of the tray
 - 10,000 alternating lateral thrusts of 150 N
 - 1000 vertical thrusts of 300 N on the tray